

# Efficiency in flow labs

Speaker:  
Detlef Claeys  
Friday 23rd, 1:00pm

## BD OneFlow™ reagents on BD FACSLyric™ Flow Cytometer: driving efficiency and cost savings

Provided since 2014 as the first ready-to-use, pre-titered, pre-dispensed, dried single-test tube format, the BD OneFlow™ reagents help in improving both laboratory processes and staff efficiencies. BD OneFlow™ reagents eliminate the need for antibody pipetting, thereby minimising operational mistakes, the risk of repeated testing and simplify antibody inventory management.

Available now also as validated assay on the BD FACSLyric™ Flow Cytometer, a new diagnostic solution combines the benefits of BD OneFlow™ reagents with the instrument sensitivity, resolution and standardisation to drive higher efficiency (reduced hands-on time and total process time) and cost savings.

[bd.com](https://www.bd.com)

BD FACSLyric™ Flow Cytometer is a Class 1 Laser Product.



BD, the BD Logo, BD FACS, OneFlow, FACSLyric and FACSuite are trademarks of Becton, Dickinson and Company or its affiliates. © 2021 BD. All rights reserved.

